

# DIGITAL Holter recorder

“SMALL”  
REDEFINED  
FOR HOLTER  
RECORDING



**H3+**

## P R O D U C T F E A T U R E S

- **Small and Lightweight** – Just 1 ounce (28 grams), so small the patient will hardly notice it is there!
- **3-Channel Recording** – Three channel recording, optional two channel bi-polar recording also available. Pacemaker spike detection standard.
- **Waveform Display** – Graphic display to preview ECG waveforms during patient hook-up.
- **Recording for 24 and 48 Hours** – With just one AAA battery the H3+ has the ability to record continuously for up to 2 days.
- **USB Data Transfer** – Recording downloaded via USB cable to H-Scribe Holter system.

## H 3 + <sup>TM</sup> D I G I T A L R E C O R D E R

Feature	Specifications
Instrument Type	Digital Holter Recorder
Input Channels	Simultaneous acquisition of two or three channels
Standard Leads Acquired	Modified I, II and V or Bipolar Channel 1 and Channel 2
Input Impedance Input Dynamic Range Electrode Offset Tolerance Frequency Response	Meet or exceed the requirements of ANSI/AAMI EC38
Digital Sampling Rate	180 s/sec/channel used for standard recording and storage.
Special Functions	Pacemaker Detection, ECG Display during hookup
A/D Conversion	12 bit
Storage	Internal, non-volatile memory
Device Classification	Type CF defibrillator proof applied parts, internally powered
Weight	1 Ounce (28 g) without battery
Dimensions	2.5 x 1.0 x .75 inches (64 x 25 x 19 mm)
Battery	1 AAA Alkaline Required for up to 48 Hours

## S T A N D A R D A C C E S S O R I E S

- Patient Cable
- Carry Case with Clip
- User Manual
- Short Form Instruction Card

### Warranty and Service

Mortara Instrument is committed to the highest level of customer support. Please contact us for the program which will be best suited for your needs.

Mortara Instrument, Inc.  
7865 North 86th Street  
Milwaukee, WI 53224  
Tel: 414.354.1600  
Tel: 800.231.7437  
Service: 888.MORTARA  
Fax: 414.354.4760

Mortara Rangoni Europe, Srl  
(European Headquarters, Italy)  
Via Cimarosa 103/105  
40033 Casalecchio di Reno (BO)  
Tel: +39.051.298.7811  
Fax: +39.051.613.3582

Mortara Instrument GmbH  
(Germany)  
Kaninenberghöhe 50  
45136 Essen  
Tel: +49.201.268311  
Fax: +49.201.268313

Mortara Instrument B.V.  
(The Netherlands)  
H. Dunantplein 6  
3731 CL De Bilt  
Postbus 131  
3720 AC Bilthoven  
Tel: +31.30.2205050  
Fax: +31.30.2201531

 **Mortara INSTRUMENT**  
The World Leader in ECG Technology

<http://www.mortara.com>  
ISO 13485 Certified