

# 12-LEAD wireless digital transmitter

THE  
WORLD'S  
ONLY  
DIAGNOSTIC,  
12-LEAD  
TRANSMITTER



## X12+

### P R O D U C T F E A T U R E S

- **Freedom of Movement** –  
The X12+ provides wireless freedom of movement for the patient and clinician during exercise stress testing.
- **Compact and Lightweight** –  
The X12+ has been designed to maximize patient comfort.
- **Diagnostic Quality** –  
Mortara wireless technology features full diagnostic ECG quality.
- **Significant Artifact Reduction** –  
The wireless solution eliminates artifact due to traditional cable movement, resulting in cleaner, higher quality tracings.
- **Comfortable ECG Acquisition for the Patient** –  
The X12+ and LeadForm patient cable are designed to eliminate the cumbersome inconvenience of heavy cabling with prepositioned leads to minimize cabling and maximize comfort.
- **The User Selectable Channels** –  
The X12+ offers a wide range of channels for interference free transmission.
- **Waveform display and Digital Menu Navigation** –  
Full graphic display to preview ECG waveforms and check lead impedance during patient hook-up, channel selection and other set up features.

## X 12 + <sup>™</sup> D I G I T A L T R A N S M I T T E R

Feature	Specifications
Instrument Type	12-lead diagnostic ECG transmitter
Input Channels	Continuous 12-lead signal acquisition and transmission
Standard Leads Acquired	I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6
Frequency Range*	608.48 to 631.52 MHz, or 2400.96 MHz to 2482.56 MHz
Special Functions	Lead impedance check, ECG display, lead fail, low battery notification, multi-purpose call, 10-wire, 5-wire and 4-wire options.
Defibrillator Protection	In compliance with AAMI standards and IEC 601-2-25
Number of Channels	Up to 256, user selectable
Function Keys	Up/right, down and enter keys for on, configuration menu navigation and call
Device Classification	Type CF defibrillator proof applied parts, battery operated
Weight	4 ounces (125g) without battery
Dimensions	2.5 x 3.5 x .98 inches (64 x 91 x 25 mm)
Battery	1 AA alkaline required

\*Operating frequency range is dependent on the X12+ part number.

## S T A N D A R D A C C E S S O R I E S

- LeadForm 12-lead Patient Cable
- Transmitter Carry Case with Belt and Strap
- User Manual
- Short Form Instruction Card
- 1 AA Alkaline Battery
- Hook-up Kit

### Warranty and Service

Mortara Instrument is committed to the highest level of customer support. Please contact us for the program which will be best suited for your needs.

Mortara Instrument, Inc.  
7865 North 86th Street  
Milwaukee, WI 53224  
Tel: 414.354.1600  
Tel: 800.231.7437  
Service: 888.MORTARA  
Fax: 414.354.4760

Mortara Rangoni Europe, Srl  
(European Headquarters, Italy)  
Via Cimarosa 103/105  
40033 Casalecchio di Reno (BO)  
Tel: +39.051.298.7811  
Fax: +39.051.613.3582

Mortara Instrument GmbH  
(Germany)  
Kaninenberghöhe 50  
45136 Essen  
Tel: +49.201.268311  
Fax: +49.201.268313

Mortara Instrument B.V.  
(The Netherlands)  
H. Dunantplein 6  
3731 CL De Bilt  
Postbus 131  
3720 AC Bilthoven  
Tel: +31.30.2205050  
Fax: +31.30.2201531

 **Mortara INSTRUMENT**  
The World Leader in ECG Technology

<http://www.mortara.com>  
ISO 13485 Certified